

# STORAGE AND DISTRIBUTION

## OVERVIEW:

- Available individually or as part of a total water treatment solution
- Meets/exceeds the exacting standards of the pharmaceutical industry and cGMP requirements
- Wide range of sizes and options available to fit your unique requirements
- Efficient component layout allows for ease of maintenance and calibration
- Customizable to customer specifications for valves, instruments and other mechanical components
- Optional complete validation package with IQ/OQ protocol execution
- Supported by MECO 24-hour customer service
- Backed by warranty

## PRODUCT HIGHLIGHTS:

### Standard Features

- Unit is completely piped, wired and assembled on skid(s) at the factory
  - Reduces installation costs and startup time
- Sanitary construction and connections
- Pumps with Variable Frequency Drives (VFD)
  - Allows for PID loop control of flow and pressure
- Hot water sanitizable
- PLC based controls

### Equipment Details

- Double tubesheet heat exchangers
- Sloped piping and low point drains
- Minimum deadlegs
- 25 Ra  $\mu\text{in}$  (0.64 Ra  $\mu\text{m}$ ) standard process contact surfaces
- Premium valves and instruments
- Polypropylene or stainless steel air lines
- Stainless steel hardware



WHERE TOMORROW GETS ITS WATER.™

