



# ES SYSTEMS by Water Works, Inc.

Ultrapure Water for Life Sciences



## 60% More Water in 50% Less Space - What are you going to do with all that extra room?

- 1 **Our design saves water & energy**  
Efficient membranes, simplified design, and maximum water recovery reduce water and power consumption.
- 2 **The integrated design combines multiple operations**  
Reduces footprint size by more than 50% and eliminates non-essential components.
- 3 **Every pipe connection is an Infrared certified weld**  
Saves inches on more than 200 connections, reducing footprint. Cleaner connections lower biofilm risk.
- 4 **Water Works customizes a UL panel to your needs**
  - PLC - Your choice
  - HMI
  - Remote Communication
- 5 **Noise Level**  
ES System - 60 dB  
Competitor – 80 dB



## ES System Benefits

The ES (Environmental Savings) Equipment line from Water Works, Inc. will exceed your expectations for ultrapure water needs.

A Water Works system design starts simply. What are your production quantity and quality requirements? Water Works unifies loop operation and make-up operation into one simple system customized to your plant.

Our systems consistently produce water that exceeds global pharmacopeia standards and can also meet Ambient Water for Injection standards. We have taken a unique approach and are a life science-focused company.

Our equipment:

- Reduces your environmental water and energy footprint increasing your company's sustainability
- Lowers your OPEX cost for chemicals, water, electricity and square footage
- Reduces the risk of biological loop contamination and potential product contamination

ES UPW Systems can handle multiple loops with flows from 5 – 100 gpm and can meet production demands up to 100 gpm.

Water Works provides a full range of on-site services, including service deionization tanks that are regenerated in a USP dedicated plant, eliminating the risk of loop contamination from co-mixing resin with an industrial or water reuse application.

**ULTRAPUREWATERWORKS.COM**  
5490 COMPLEX STREET, STE #602  
SAN DIEGO, CA 92123  
(858) 499-0119

